



TA 20 Amplified Telecoils Coding Scheme

The purpose of the color coding process is to improve product identification for all amplified telecoils and we have therefore added two color stripes on the hybrid encapsulate.

Description of Coding

Since Sonion offers a variety of amplified telecoils with potentially the same physical appearance, we have implemented a procedure to add two color stripes to the top of the black encapsulate. Currently, we have identified the use of six different colors that will best differentiate the telecoils.

Each amplified telecoil part number has been assigned a different color code. To protect the color marking from physical wear, we implement the color marking procedure before the top layer protective coating is applied.

We are providing this service to only our amplified telecoils and not to our passive coils. Passive coils have a variety of physical and electrical differences to help you more easily identify the part. As new amplified telecoil models are developed, we will continue to create new color code schemes to help you better identify the product.

Technical Bulletin

Results

For the complete understanding to the color coding matrix, please see the below overview.

Amplified Telecoil Color Code Matrix

Item no.	Type Name	Color Marking Dot 1	Dot 2
9010-2000180	TA 20 AD 01	yellow	blue
9010-2003066	TA 20 AD 03	white	blue
9009-3005110	TA 20 AD 11	yellow	no line
9008-3009565	TA 20 AD 21	blue	blue
9010-2000341	TA 20 AE 01	yellow	yellow
9010-2002827	TA 20 AE 03	yellow	green
9010-3006102	TA 20 AE 10	white	red
9009-3005108	TA 20 AE 11	yellow	red
9008-3008312	TA 20 AE 21	blue	grey
9010-3002730	TA 20 BE 01	green	white
9008-3009894	TA 20 BE 21	blue	yellow
9010-3002733	TA 20 CE 01	green	grey
9010-3010785	TA 20 CE 03	white	grey
9008-3008328	TA 20 CE 21	blue	red
9010-3002276	TA 20 DE 01	yellow	white
9008-3016674	TA 32 CE 21	white	white

Item no.	Type Name	Dot 1	Dot 2
9010-2000180	TA 20 AD 01	yellow	blue



Item no.	Type Name	Dot 1	Dot 2
9010-2003066	TA 20 AD 03	white	blue



Item no.	Type Name	Dot 1	Dot 2
9009-3005110	TA 20 AD 11	yellow	no line



Technical Bulletin

Item no.

9008-3009565

Type Name

TA 20 AD 21

Dot 1

blue

Dot 2

blue

**Item no.**

9010-2000341

Type Name

TA 20 AE 01

Dot 1

yellow

Dot 2

yellow

**Item no.**

9010-2002827

Type Name

TA 20 AE 03

Dot 1

yellow

Dot 2

green

**Item no.**

9010-3006102

Type Name

TA 20 AE 10

Dot 1

white

Dot 2

red

**Item no.**

9010-3005108

Type Name

TA 20 AE 11

Dot 1

yellow

Dot 2

red

**Item no.**

9008-3008312

Type Name

TA 20 AE 21

Dot 1

blue

Dot 2

grey



Technical Bulletin

Item no.

9010-3002730

Type Name

TA 20 BE 01

Dot 1

green

Dot 2

white

**Item no.**

9008-3009894

Type Name

TA 20 BE 21

Dot 1

blue

Dot 2

yellow

**Item no.**

9010-3002733

Type Name

TA 20 CE 01

Dot 1

green

Dot 2

grey

**Item no.**

9010-3010785

Type Name

TA 20 CE 03

Dot 1

white

Dot 2

grey

**Item no.**

9008-3008328

Type Name

TA 20 CE 21

Dot 1

blue

Dot 2

red

**Item no.**

9010-3002276

Type Name

TA 20 DE 01

Dot 1

yellow

Dot 2

white



Technical Bulletin

Item no.

9008-3016674

Type Name

TA 32 CE 21

Dot 1

white

Dot 2

white

**Remark**

Please contact your local sales representative for any questions you may have. We thank you for your business and are pleased to continually improve the quality and performance of our products.